Pictogrammes/Pictogrammes/Pictogrammen :

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ENGLISH- CYRANOSE® GLOBAL SYSTEM INSTRUCTIONS FOR USE

CEREDAS, a French designer, manufacturer and distributor of medical devices since 1993, and certified as per ISO standard 13485:2016, has a range of medical devices: CYRANOSE® Global System dedicated to the comfort of patients with laryngectomy and tracheotomy, which includes patented products: the CYRANOSE® artificial nose or HME and the TRACHEOCLEAN® disposable tracheal protection, both reimbursed by the national healthcare system in France and Belgium.

CEREDAS expressly invites you to carefully read the instruction given in this leaflet, as well as the labelling of these same Medical Devices.

The CEREDAS medical device "CYRANOSE® HME*" from the "CYRANOSE® Global System™ range" is a class I medical device as per CE marking - Regulation: MDR EU 2017/745 -, which underwent several technical studies on heating, moistening and filtration of inhaled air by the National Testing Laboratory (LNE) prior to its Marketing Authorisation - ref. Technical CE Marking File -; The CYGRIPONE™ and CYTUBEMAJOR™ laryngectomy tubes and tracheostoma supports are class IIb devices.

The technical and medical data that characterise the performance of the medical device as mentioned above and claimed before the Competent Authorities, have been the subject of several technical and medical publications in the history of the medical device in question, and in particular by the French National Authority for Health (HAS) as part of the validation by the CNEDIMTS of the "list of reimbursed products and service referencing in France".

If a serious incident is identified, it must be reported to the Manufacturer CEREDAS (ceredas@ceredas.com) and the Competent Authority (ANSM in France: <https://ansm.sante.fr/>)

“serious incident”, any incident that directly or indirectly caused, could have caused or could cause:

- a) the death of a patient, of a user or of any other person;***
- b) serious deterioration, temporary or permanent, in the health condition of a patient, of a user or of any other person;***
- c) a serious threat to public health)***

USE OF THE DEVICES

CLAIMED USE

The CYTUBEMAJOR™ and CYGRIPONE™ devices are intended to be inserted in the tracheostoma of patients with laryngectomy in order to calibrate and maintain the opening.

The CYRANOSE® device is intended to moisten, warm and filter the air inhaled by the patient and makes it possible to increase pulmonary resistance. It can be placed directly on the patient's skin (with or without the CYGRIPONE™) or fixed on the base of the CYTUBEMAJOR™ device.

CLINICAL CLAIM

The clinical performance, making it possible to produce the "Clinical Benefit", is related to the indirect effect of the technical characteristics of the CYRANOSE®, CYTUBEMAJOR™ and CYGRIPONE™ devices. According to the principle of how the device works, the technical characteristics are as follows: CYTUBEMAJOR™ and CYGRIPONE™

- Ability to maintain the calibre of the tracheostoma due to - mechanically suitable and biocompatible materials;
 - A diameter and length adapted to the patient's tracheostoma.
- CYTUBEMAJOR™
- Ability to act as a support for the CYRANOSE® device via its structure.
- CYRANOSE®
- Ability to moisten, heat and filter inhaled air due to its structure and the materials of the various parts that it is composed of.

CLINICAL BENEFIT

CYTUBEMAJOR™ and CYGRIPONE™ laryngectomy tubes enable the calibration and protection of the tracheostoma of patients with laryngectomy due to their form and their composition. They can also be used alone or in combination with the CYRANOSE® artificial nose or HME (Heat and Moisture Exchanger) and thus contribute to the same clinical benefit as the CYRANOSE® device. CYRANOSE® is a reusable medical device intended to heat, moisten and filter inhaled air in order to improve the breathing of patients with laryngectomy. In all seasons and in all places, the person with a laryngectomy experiences first hand the changes in the surrounding air: dryness in the summer but also in the winter with certain heating systems, sudden cooling in winter but also in the summer in air conditioned rooms, pollen in the spring, etc. As the inhaled air is not filtered, heated or moistened by the upper respiratory tract, the respiratory system is particularly exposed and can have the following consequences:

- hypersensitivity of the epithelium,
- pollution and inhalation of particles,
- cough,
- increased volume of mucus,
- dryness,
- irritation, mucus plugs and crusts.

CYRANOSE®, used alone or on the CYTUBEMAJOR™ and CYGRIPONE™ devices, re-establishes a conditioning of the air that is close to the natural system of breathing. It enables:

- filtration of the air with a performance equivalent to that of the upper airways;
- a reduction in cough and therefore secretions;
- moistening of the airways, preventing the formation of mucus crusts and plugs.
- improved breathing.

55.3.1 CONTRAINDICATIONS

The contraindications of the CYRANOSE heat and moisture exchanger (HME) are as follows:

- In case of intolerance of the double-sided adhesive;
 - When the morphology of the tracheostoma does not enable the proper positioning of the base of the "CYRANOSE heat and moisture exchanger (HME)" respiratory prosthesis (leaks).
- In these 2 cases, the patient can use the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support

55.3.2 INDICATIONS

Tracheostoma device for patients with or without a voice prosthesis after total laryngectomy or total laryngopharyngectomy.

55.3.3 TARGET POPULATION

Patients with laryngectomy or laryngopharyngectomy with or without a voice prosthesis.

55.3.4 USERS

Healthcare professionals in hospital settings and patients with laryngectomy or laryngopharyngectomy.

CYRANOSE ECH


8.4.5.4.2.1 • RECOMMENDATIONS AND REMINDERS BEFORE PLACING THE CYRANOSE® ARTIFICIAL NOSE OR HME*

- First carefully wash your hands
- Do not apply the adhesive to irritated, sensitive, fragile or unhealed skin

For persons with a risk of skin reactions, the PS (sensitive skin) adhesive is recommended after testing

- Always clean your skin before applying a new adhesive

Strict hygiene is necessary for the proper maintenance of the device: use a wash cloth to clean the skin around the tracheostoma with warm water and pH neutral soap before applying the adhesive. Wipes exist for persons subject to frequent detachment to help improve adherence.

PLACEMENT OF THE CYRANOSE® ARTIFICIAL NOSE OR HME WITH ITS SUPPORT, BASE AND ADHESIVE.

STEP 1: SEPARATION OF THE ADHESIVE FROM ITS PROTECTION

Detach the first protective layer from the adhesive (A) (split part)

STEP 2: PLACEMENT OF THE ADHESIVE ON THE ADHESIVE GUIDE

Place the adhesive, with the sticky side towards you, on the blue adhesive guide

STEP 3: PLACEMENT OF THE BASE ON THE ADHESIVE

Position and stick the base (B) on the adhesive using the blue adhesive guide. Remove it while continuing to firmly apply the adhesive on the base by pressing with the thumb and index finger.

STEP 4: REMOVAL OF THE SECOND PROTECTION

Detach the second protection using the tab located on the side of the adhesive (facilitate this by first removing the rest of the first protection and the adhesive located on the tab).

STEP 5: PLACEMENT OF THE BASE AROUND THE TRACHEOSTOMA

Stick the base (B) around the tracheostoma by gently pulling on the skin and using the position that you find most suitable. Properly apply the adhesive on the base and around the tracheostoma using your thumb to massage the parts for 10 seconds in order to make sure that the adhesive fully adheres to the skin.

STEP 6: PLACEMENT OF THE SECRETION TRAP IN THE BASE


Place the secretion trap (C) in the base while holding the picot edge between the thumb and index finger. The secretion trap should not touch the edge of the tracheostoma or any part due to a risk of irritation.

STEP 7: PLACEMENT OF THE FOAM FILTER IN THE CASING

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STEP 8: CLOSURE OF THE DEVICE

Close the device by placing the casing (E) on the base ring (B), turn to the left until it enters the notches and then tighten gently in order to block it.

MANAGEMENT OF MUCUS WHEN COUGHING

Always hold the device in place so as not to weaken the quality of the adhesion. For use in a hospital: the mucus can be removed from the secretion trap by suction. For ambulatory use: always carry a secretion trap and a foam filter with you so that you can change the ones soiled with mucus.

IN ASSOCIATION WITH A VOICE PROSTHESIS

If the stopper housing is used with a voice prosthesis, the placement of the base and the adhesive is important in order to avoid air leaks under the adhesive. Check that the adhesive sticks well around the entire tracheostoma. To make sure that the stopper housing is well closed when you use the voice prosthesis, place two or three fingers on the cover to ensure a good seal around the silicone joint during phonation.

Due to the expiratory pressure used to operate the valve, after a few hours a leak may occur around the tracheostoma, preventing proper operation of the voice prosthesis and the valve. Put the stopper housing back in place and insert the automatic valve or apply a fresh adhesive.

IN ASSOCIATION WITH THE CYADAPTOR™ ADAPTOR AND THE CYVALVEFREEDOM™ HANDS-FREE VALVE

Due to the expiratory pressure used to operate the valve, after a few hours a leak may occur around the tracheostoma, preventing proper operation of the voice prosthesis and the valve. Put the stopper housing back in place and insert the automatic valve or apply a fresh adhesive.

IN ASSOCIATION WITH THE TRACHEODOUCHE™ SHOWER PROTECTOR

If the TRACHEODOUCHE™ shower protector is used, remove all of the elements from the prosthesis and only keep the base attached to the skin with the adhesive to attach the TRACHEODOUCHE™ with the opening towards the bottom.

N.B.: when taking a shower, check for proper sealing between the adhesive and the skin.

IN ASSOCIATION WITH THE SILICONE CYGRIPONE™ LARYNGECTOMY TUBE

The flange of the CYGRIPONE™ laryngectomy tube is compatible with the diameter of the base of the CYRANOSE® artificial nose or HME. If the CYGRIPONE™ laryngectomy tube is used with the CYRANOSE® artificial nose or HME, the CYRANOSE® base and adhesive must be positioned (centred) beforehand around the tracheostoma. The CYGRIPONE™ laryngectomy tube must then be inserted in the tracheostoma. In case of cough, if the entire device is congested, remove the housing, the secretion trap, clean the mucus in the base with a tissue, then use the grip tab to remove the CYGRIPONE™ laryngectomy tube and clean it.

IN ASSOCIATION WITH THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

The CYRANOSE® artificial nose or HME can be worn using the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support after total laryngectomy or when the patient needs to allow his or her skin to rest.

PLACEMENT OF THE CYRANOSE® ARTIFICIAL NOSE OR HME ON THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

Carefully wash your hands beforehand, then take the following from the Cyranose® starter kit: 1 secretion trap (C) + 1 foam filter (D) + 1 housing (E), which have been previously disassembled from the base, with everything packaged in a translucent bag in the starter kit.

STEP 1: PLACEMENT OF THE SECRETION TRAP

Place the secretion trap (C) in the ring of the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support while holding the picot edge between the thumb and index finger.

STEP 2: PLACEMENT OF THE FOAM FILTER (D) IN THE HOUSING (E)

STEP 3: CLOSURE OF THE DEVICE

Close the device by placing the housing (E) on the ring of the CYTUBEMAJOR laryngectomy tube and tracheostoma support, turn to the left until it enters the notches and then tighten gently in order to block it.

REMINDERS REGARDING THE USE OF THE CYRANOSE® ARTIFICIAL NOSE OR HME

The morphology of the neck can change over time, so it is necessary to select the base model that is best adapted to this morphology. If no adhesive support is possible or satisfactory, the device may be worn with the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support. The quality of the adherence is not guaranteed with heavy sweating. With sustained effort and/or when the outside temperature is very low, the difference in temperature may cause strong condensation on the metal parts and induce detachment of the device. With sustained effort and/or when the outside temperature is very low, the difference in temperature may cause strong condensation on the metal parts and saturate the foam filter with water. It may then be sponged up or changed with a clean and dry filter. It is necessary to check each element every day, in particular the metal parts of the device (traces of oxidation or corrosion).

REMOVAL OF THE CYRANOSE® ARTIFICIAL NOSE OR HME WITH ITS SUPPORT, BASE AND ADHESIVE FIG. 5

The adhesive should not remain on the skin for more than 24 hours. Allow some time between each removal of the device.

First carefully scrub your hands.

- 1- Remove the elements from the device: housing, filter, secretion trap
- 2- Clean the edges of the tracheostoma if there is any mucus using a clean tissue
- 3- Take one side of the base and pull gently by rolling the adhesive until the unit is slowly removed (adhesive + base). Remove the base from the skin by pulling gently.
- 4- Remove the adhesive using the pre-cut notch inside the base and discard in the garbage bin

CLEANING THE CYRANOSE® ARTIFICIAL NOSE OR HME

The elements in the device should be cleaned after each use. The adhesive should be removed from its support after each use and discarded in the garbage bin. Base (B): Clean with warm soapy water, then rinse and dry in the folds of a clean and dry towel. Do not use detergents or chemical products or alcohol as this may destroy the adhesive properties of the base on the adhesive. Secretion trap (C): If the mucus has dried on the mesh and is too hard to remove, soak the secretion trap in clean water for several hours and then clean with warm soapy water with a soft brush. Rinse and dry in the folds of a clean and dry towel. Foam filter (D): Clean with warm soapy water and rinse thoroughly and dry in the folds of a clean and dry towel. The foam filter may be used after each washing for 30 days. After this time, discard it in the garbage bin. The filter may also be used for single use (no cleaning). Standard housing (E): Clean the housing with warm soapy water with a soft brush, especially the notched inner part and/or using a compress with 70% alcohol. Rinse and dry in the folds of a clean and dry towel. Stopper housing (E): As with the standard housing, being careful to remove the stainless steel spring and not damage the seal during cleaning with a soft brush. Rinse and dry in the folds of a clean and dry towel.

STORAGE

The Cyranose artificial nose or HME should be stored protected from direct light in a room with a temperature between: 10–27°C / 50–80°F and humidity between 40–60%.

RESIDUAL RISK
In case of allergic reactions (redness, itching, irritation), it is recommended to check whether cleaning has been performed as per the recommendations in section “Cleaning the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support” (or section “Cleaning the CYGRIPONE™ laryngectomy tube”). If this is the case, you should stop using the CYTUBEMAJOR™ (or CYGRIPONE™) devices, consult the doctor who prescribed these medical devices and inform CEREDAS by phone, letter or email.

REMOVAL OF THE DOUBLE-SIDED ADHESIVE

The double-sided adhesive is to be used only once. The double-sided adhesive should be removed from its support after each use and discarded in the garbage bin.

CYRANOSE® ARTIFICIAL NOSE OR HME - STANDARD STARTER KIT including: 1 standard housing, 3 secretion traps, 2 bases, 6 filters, 20 adhesives, 1 adhesive guide				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
A, standard	CYR1013	2447403	100%	CYRANOSE HME 1 starter kit
A, sensitive skin	CYR1014			
D, standard	CYR1017			
D, sensitive skin	CYR1018			
G, standard	CYR1019			
G, sensitive skin	CYR1020			
E, standard	CYR1021	2430331	100%	CYRANOSE HME 1 starter kit
E, sensitive skin	CYR1022			
CYRANOSE® ARTIFICIAL NOSE OR HME - STOPPER STARTER KIT including: 1 stopper housing, 3 secretion traps, 2 bases, 6 filters, 20 adhesives, 1 adhesive guide				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
A, standard	CYR2013	2447834	100%	CYRANOSE HME 1 starter kit
A, sensitive skin	CYR2014			
D, standard	CYR2017			
D, sensitive skin	CYR2018			
G, standard	CYR2019			
G, sensitive skin	CYR2020			
E, standard	CYR2021	2416242	100%	CYRANOSE HME 1 starter kit
E, sensitive skin	CYR2022			
CYRANOSE® ARTIFICIAL NOSE OR HME - REPLACEMENT PARTS Replacement kit including: 100 adhesives and 3 filters				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
A, standard	CYR3001	2482960	100%	2 replacement kits
A, sensitive skin	CYR3002			
D, standard	CYR3005			
D, sensitive skin	CYR3006			
G, standard	CYR3007			
G, sensitive skin	CYR3008			
E, standard	CYR3009	2487549	100%	2 replacement kits
E, sensitive skin	CYR3010			
CYRANOSE® ARTIFICIAL NOSE OR HME - ADHESIVE REPLACEMENT PARTS (Pack of 100 adhesives)				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
A, standard	CYR4029	2450894	100%	CYRANOSE HME 200 adhesives + Form + Standard or Sensitive Skin
A, sensitive skin	CYR4030			
D, standard	CYR4031			
D, sensitive skin	CYR4032			
G, standard	CYR4033			
G, sensitive skin	CYR4034			
E, standard	CYR4035	2475902	100%	CYRANOSE HME 200 adhesives + Form + Standard or Sensitive Skin
E, sensitive skin	CYR4036			

CYRANOSE® ARTIFICIAL NOSE OR HME - REPLACEMENT PARTS FOAM FILTER (Pack of 3 filters)				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
	CYR4016	2492450	100% 3 filters every 3 months	CYRANOSE HME 3 filters
CYRANOSE® ARTIFICIAL NOSE OR HME - BASE REPLACEMENT PARTS				
Base, Adhesive	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
A	CYR4009	2469273	100% 1 per year	CYRANOSE HME 1 Reinforced base + Form
D	CYR4011			
G	CYR4014			
E	CYR4012			
CYRANOSE® ARTIFICIAL NOSE OR HME - SECRETION TRAP REPLACEMENT PARTS				
	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
	CYR4015	2401654	100% 1 per year	CYRANOSE HME 1 secretion trap
CYRANOSE® ARTIFICIAL NOSE OR HME - CASING REPLACEMENT PARTS				
	Réf.	LPPR [List of reimbursed products and services] code	Reimbursement by national health insurance system	Name to be entered
Standard	CYR4001	2460438	100% 1 every 5 years	CYRANOSE HME 1 standard housing
Stopper	CYR4002	2419737	100% 1 every 5 years	CYRANOSE HME 1 stopper housing

TUBE AND TRACHEOSTOMA SUPPORT CYTUBEMAJOR

PLACEMENT OF THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

STEP 1: PRIOR RINSING OF THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

Carefully wash your hands before rinsing the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support

STEP 2: PLACEMENT OF THE CANNULA ATTACHMENT IN THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT.

After having selected the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support model (diameter, length, fenestrated or non-fenestrated), place the cannula attachment in the side openings of the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support

STEP 3: INSERTION OF THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT IN THE TRACHEOSTOMA

Gently insert the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support in the tracheostoma, while holding the flange on both sides between the thumb and index finger, until it reaches the edges of the tracheostoma.

To facilitate placement in the tracheostoma, depending on your prescribing physician's instructions, use silicone oil. Do not use oil to facilitate insertion or the silicone could be damaged.

If the patient has a voice prosthesis, the patient or healthcare professional must take care not to hit it when inserting the device into the tracheostoma.

STEP 4 CLOSURE OF THE CANNULA ATTACHMENT

Tighten the cannula attachment slightly around the neck in order to hold the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support in place

THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT IN ASSOCIATION WITH A TRACHEAL COMPRESS

A tracheal compress can be placed between the patient's skin (depending on the condition of the skin) and the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support. Insert the laryngectomy tube into the perforated part in the middle of the compress for this purpose. Then insert the laryngeal tube into the patient's tracheostoma.

REMOVAL OF THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

After carefully washing your hands, proceed with the removal of the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support by holding each side of the flange between the thumb and index finger, after having detached the cannula attachment beforehand.

CLEANING OF THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

First carefully scrub your hands. Then clean the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support with a soft brush and soapy water and then rinse thoroughly. Use a pH neutral soap. Important: do not use detergents or chemical products.

The colour of the CYTUBEMAJOR™ may change over time and turn slightly yellow. To facilitate placement in the tracheostoma, depending on your prescribing physician's instructions, use silicone oil. Non-sterile device.

Store the CYTUBEMAJOR™ in a dark, clean and dry place.

RESIDUAL RISK


In case of allergic reactions (redness, itching, irritation), it is recommended to check whether cleaning has been performed as per the recommendations in section "Cleaning the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support" (or section "Cleaning the CYGRIPONE™ laryngectomy tube"). If this is the case, you should stop using the CYTUBEMAJOR™ (or CYGRIPONE™) devices, consult the doctor who prescribed these medical devices and inform CEREDAS by phone, letter or email.

CONTRAINDICATION

- Use of any form of mechanical ventilation
- Tissue damage in the trachea or tracheostoma;

ELIMINATION OF CYTUBEMAJOR™ WASTE

At the end of its life (6 months), this material should be disposed of with the household waste in accordance with your local waste sorting instructions.

CYTUBEMAJOR™: LARYNGECTOMY TUBE AND Tracheostoma SUPPORT - STANDARD				
Size 8 int. Ø 9.5 mm ext. Ø 12 mm	Ref.	LPPR code	Reimbursement by national health insur- ance system	Name to be entered
L 27mm	CBC0827	2471212	100 % 2 per year	2 CYTUBEMAJOR™ standard + size + length or reference
L 36mm	CBC0836			
L 55mm	CBC0855			
Taille 9 Øint 10,5mm Øext 13,5mm	Réf.			
L 27mm	CBC0927			
L 36mm	CBC0936			
L 55mm	CBC0955			
Taille 10 Øint 12mm Øext 15mm	Réf.			
L 27mm	CBC1027			
L 36mm	CBC1036			
L 55mm	CBC1055			
Taille 12 Øint 13,5mm Øext 17mm	Réf.			
L 27mm	CBC1227			
L 36mm	CBC1236			
L 55mm	CBC1255			

CYTUBEMAJOR™: LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT - FENESTRATED				
Size 8 int. Ø 9.5 mm ext. Ø 12 mm	Ref.	LPPR code	Reimbursement by national health insur- ance system	Name to be entered
L 27mm	CBCF0827	2407728	100 % 2 per year	2 CYTUBEMAJOR™ fenestrated + size + length or reference
L 36mm	CBCF0836			
L 55mm	CBCF0855			
Taille 9 Øint 10,5mm Øext 13,5mm	Réf.			
L 27mm	CBCF0927			
L 36mm	CBCF0936			
L 55mm	CBCF0955			
Taille 10 Øint 12mm Øext 15mm	Réf.			
L 27mm	CBCF1027			
L 36mm	CBCF1036			
L 55mm	CBCF1055			
Taille 12 Øint 13,5mm Øext 17mm	Réf.			
L 27mm	CBCF1227			
L 36mm	CBCF1236			
L 55mm	CBCF1255			

PLACEMENT OF THE CYGRIPONE™ LARYNGECTOMY TUBE

INSERTION RECOMMENDATIONS: The colour of the CYGRIPONE™ may change over time and yellow. To facilitate placement in the tracheostoma, depending on your prescribing physician's instructions, use silicone oil. Non-sterile device.

STEP 1: PRIOR RINSING OF THE TUBE

Rinse the CYGRIPONE™ laryngectomy tube beforehand. The CYGRIPONE™ laryngectomy tube is put in place by the healthcare professional or the patient him or herself under the supervision of a healthcare professional.

STEP 2: PLACEMENT WITH OR WITHOUT THE CANNULA ATTACHMENT

The CYGRIPONE™ laryngectomy tube can be put in place:
 - With a cannula attachment that must be inserted in the side openings of the laryngeal tube.
 - Without a cannula attachment and then held in place by the CYRANOSE® artificial nose or HME with its base and adhesive

STEP 3: INSERTION OF THE TUBE IN THE TRACHEOSTOMA

Gently insert the CYGRIPONE™ laryngectomy tube in the tracheostoma, while holding the flange on both sides between the thumb and index finger, until it reaches the edges of the tracheostoma. If the patient has a voice prosthesis, the healthcare professional must take care not to hit it when inserting the device into the tracheostoma.

REMOVAL OF THE CYGRIPONE™ LARYNGECTOMY TUBE

Pull gently using the pull tab located above the flange and proceed with the removal of the CYGRIPONE™ laryngectomy tube from the tracheostoma.

CLEANING OF THE CYGRIPONE™ LARYNGECTOMY TUBE

First carefully scrub your hands. Then clean the CYGRIPONE™ with a soft brush and soapy water and thoroughly rinse. Use a pH neutral soap. Important: do not use detergents or chemical products.



STORAGE
Store the CYGRIPONE™ in a dark, clean and dry place.

RESIDUAL RISK
In case of allergic reactions (redness, itching, irritation), it is recommended to check whether cleaning has been performed as per the recommendations in section "Cleaning the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support" (or section "Cleaning the CYGRIPONE™ laryngectomy tube"). If this is the case, you should stop using the CYTUBEMAJOR™ (or CYGRIPONE™) devices, consult the doctor who prescribed these medical devices and inform CEREDAS by phone, letter or email.

ELIMINATION OF CYGRIPONE™ WASTE

At the end of its life (6 months), this material should be disposed of with the household waste in accordance with your local waste sorting instructions.


CYGRIPONET™ LARYNGECTOMY TUBE STANDARD				
Size 8 int. Ø 9.5 mm ext. Ø 12 mm	Ref.	LPPR code	Reimbursement by national health insur- ance system	Name to be entered
L 18mm	CNSF0818	2436606	100% 2 per year	2 CYGRIPONET™ STANDARD + SIZE + LENGTH OR REFERENCE
L 27mm	CNSF0827			
L 36mm	CNSF0836			
L 55mm	CNSF0855			
Taille 9 Øint 10,5mm Øext 13,5mm	Réf.			
L 18mm	CNSF0918			
L 27mm	CNSF0927			
L 36mm	CNSF0936			
L 55mm	CNSF0955			
Taille 10 Øint 12mm Øext 15mm	Réf.			
L 18mm	CNSF1018			
L 27mm	CNSF1027			
L 36mm	CNSF1036			
L 55mm	CNSF1055			
Taille 12 Øint 13,5mm Øext 17mm	Réf.			
L 18mm	CNSF1218			
L 27mm	CNSF1227			
L 36mm	CNSF1236			
L 55mm	CNSF1255			

CYGRIPONET™ FENESTRATED LARYNGECTOMY TUBE				
Size 8 int. Ø 9.5 mm ext. Ø 12 mm	Ref.	LPPR code	Reimbursement by national health insur- ance system	Name to be entered
L 18mm	CNSF0818	2498748	100% 2 per year	2 CYGRIPONET™ fenestrated + size + length or reference
L 27mm	CNSF0827			
L 36mm	CNSF0836			
L 55mm	CNSF0855			
Taille 9 Øint 10,5mm Øext 13,5mm	Réf.			
L 18mm	CNSF0918			
L 27mm	CNSF0927			
L 36mm	CNSF0936			
L 55mm	CNSF0955			
Taille 10 Øint 12mm Øext 15mm	Réf.			
L 18mm	CNSF1018			
L 27mm	CNSF1027			
L 36mm	CNSF1036			
L 55mm	CNSF1055			
Taille 12 Øint 13,5mm Øext 17mm	Réf.			
L 18mm	CNSF1218			
L 27mm	CNSF1227			
L 36mm	CNSF1236			
L 55mm	CNSF1255			

OTHER OPTIONAL COMPONENTS CAN ALSO BE COMBINED WITH CYRANOSE HME

ADAPTATOR CYADAPTOR™

PLACEMENT OF THE CYADAPTOR™ ADAPTOR

After having put the foam filter of the CYRANOSE® artificial nose or HME inside the CYADAPTOR™ adaptor, place the adaptor on the ring of the base of the CYRANOSE® artificial nose or HME, then turn to the left until it enters the notches and then tighten gently in order to block it.

REMOVAL OF THE CYADAPTOR™ ADAPTOR

Turn the CYADAPTOR™ universal adaptor to the right to the stop of the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support or the base of the CYRANOSE® artificial nose or HME, then remove it from the notches.

CLEANING THE CYADAPTOR™ UNIVERSAL

ADAPTOR

First carefully scrub your hands. Then clean the CYADAPTOR™ with a soft brush and soapy water and thoroughly rinse. Use a pH neutral soap. Important: do not use detergents or chemical products.


STORAGE

The CYADAPTOR™ adaptor should be stored in a dark, clean and dry place.

ELIMINATION OF CYADAPTOR™ WASTE

At the end of its life (6 months), this material should be disposed of with the household waste in accordance with your local waste sorting instructions.

HANDS-FREE VALVE CYVALVEFREEDOM™

PLACEMENT OF THE CYVALVEFREEDOM™ HANDS-FREE VALVE (VIA THE BASE AND THE ADHESIVE OF THE CYRANOSE HME)

Generally a few months (6 to 9 months) after the end of radiotherapy and depending on the condition of the patient's skin as well as the morphology around the tracheostoma, the use of the CYVALVEFREEDOM™ hands-free valve will be indicated in combination with the components of the CYRANOSE® artificial nose or HME: Adhesive, base, secretion trap and foam filter. In order to work correctly, the proper seal around the adhesive of the CYRANOSE® artificial nose or HME must be checked. The CYVALVEFREEDOM™ hands-free valve is thus maintained on the base of the CYRANOSE® artificial nose or HME (and by means of the adhesive).

USE OF THE CYVALVEFREEDOM™ HANDS-FREE VALVE

Due to the expiratory pressure used in order to operate the valve, after a few hours a leak may occur around the tracheostoma, preventing the proper operation of the implant and the valve. It is then possible to temporarily use a finger to plug the tracheostoma before changing: - either the adhesive support of the CYRANOSE® artificial nose or HME - or the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support. When coughing, the silicone membrane comes out of the valve in order to let the air pass. It must then be put back in place inside the valve.

ADJUSTING THE CYVALVEFREEDOM™ HANDS-FREE VALVE

Adjust the opening of the CYVALVEFREEDOM™ hands-free valve by turning it to the right or to the left. The adjustments are different for each person and must be carried out on yourself. For a maximum opening, the valve should be turned to your right and for a minimum opening, it must be turned to your left. The smaller the opening, the more sensitive the closure of the membrane will be.

IN ASSOCIATION WITH THE CYRANOSE® ARTIFICIAL NOSE OR HME

The CYVALVEFREEDOM™ hands-free valve is positioned on the base and adhesive of the CYRANOSE® artificial nose or HME. In combination with the secretion trap and foam filter, the device enables the user to speak without using his or her hands and benefit from the advantages of the CYRANOSE® artificial nose or HME.

REMOVAL OF THE CYVALVEFREEDOM™ HANDS-FREE VALVE

The CYVALVEFREEDOM™ hands-free valve is removed: Either by gently pulling outwards on its CYADAPTOR™ support while maintaining the adaptor with one hand, Or by removing the CYADAPTOR™+CYVALVEFREEDOM™ together by turning the CYADAPTOR™ to the patient's right in order to remove the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support.

CLEANING OF THE CYVALVEFREEDOM™ HANDS-FREE VALVE

First carefully scrub your hands. Then clean the CYVALVEFREEDOM™ hands-free valve with a soft brush and soapy water and then thoroughly rinse. Use a pH neutral soap. Detergents and chemical products must absolutely not be used.

The CYVALVEFREEDOM™ hands-free valve should be stored in a dark, clean and dry place.

ELIMINATION OF CYVALVEFREEDOM™ WASTE

At the end of its life (1 year), this material should be disposed of with the household waste in accordance with your local waste sorting instructions.

HANDS-FREE KIT: CYADAPTOR™ (universal adaptor) + CYVALVEFREEDOM™ (hands-free valve)				
Name	Ref.	LPPR code	Reimbursement by national health insurance system	Name to be entered
Universal adaptor	AUD4001	2420806	100% 2 per year	1 CYVALVEFREEDOM + 2 CYADAPTOR
Hands-free valve	KML4001	2449201	100% 1 per year	

LE TRACHEODOUCHE™ SHOWER PROTECTOR

GENERAL DESCRIPTION

The TRACHEODOUCHE™ shower protector is part of the CYRANOSE® GLOBAL SYSTEM according to the decision tree

USE OF THE TRACHEODOUCHE™

The TRACHEODOUCHE™ shower protector can be used depending on the various stages in the care process of the patient with laryngectomy:

-After the operation: When healing is considered sufficient, the patient can begin taking showers, under the supervision of nursing personnel, and using the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support, and by checking that the seal between the tracheostoma and the CYTUBEMAJOR™ is satisfactory

-On a daily basis: When it is no longer necessary to use the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support, TRACHEODOUCHE™ can be used with the adhesive and base of the CYRANOSE® artificial nose or HME.

PLACEMENT OF THE TRACHEODOUCHE™ SHOWER PROTECTOR

The shower protector can be used with the adhesive (A) and the base (B) of the CYRANOSE® artificial nose or HME. Remove all of the parts of the CYRANOSE® artificial nose or HME and only keep the base attached with the adhesive on the skin, then insert the TRACHEODOUCHE™ shower protector with the opening towards the bottom in the ring of the CYRANOSE® artificial nose or HME. Do not forget to check the seal between the adhesive and the skin.

PLACEMENT OF THE TRACHEODOUCHE™ WITH THE CYTUBEMAJOR™ LARYNGECTOMY TUBE AND TRACHEOSTOMA SUPPORT

The TRACHEODOUCHE™ shower protector is generally used about 10 days after the laryngectomy when the skin has healed. The TRACHEODOUCHE™ shower protector is placed on the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support by the healthcare professional or by the patient him or herself under the supervision of the healthcare professional and as per the care protocols used in the hospital department. The TRACHEODOUCHE™ shower protector must be inserted in the ring of the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support. It is necessary to check that the seal between the tracheostoma and the CYTUBEMAJOR™ laryngectomy tube and tracheostoma support is satisfactory

REMOVAL OF THE TRACHEODOUCHE™ SHOWER PROTECTOR

The TRACHEODOUCHE™ shower protector is removed by pulling it gently outwards.

CLEANING OF THE TRACHEODOUCHE™ SHOWER PROTECTOR

First carefully scrub your hands. Then clean the TRACHEODOUCHE™ shower protector with a soft brush and soapy water and then thoroughly rinse and dry. Use a pH neutral soap. Important: do not use detergents or chemical products.

The TRACHEODOUCHE™ shower protector should be stored in a dark, clean and dry place.

ELIMINATION OF TRACHEODOUCHE™ WASTE

At the end of its life (1 year), this material should be disposed of with the household waste in accordance with your local waste sorting instructions.

TRACHEODOUCHE™ - SHOWER PROTECTOR				
	Ref.	LPPR code	Reimbursement by national health insurance system	Name to be entered
	TRD4001	2444876	100% 1 per year	1 TRACHEODOUCHE

TRACHEOCLEAN

DESCRIPTION - The TRACHEOCLEAN® disposable tracheal protection is a disposable tracheal protection made of non-woven material intended for persons with laryngectomy and persons with tracheotomy. It is attached around the neck with two cords (1), one of which is pre-cut (2) in order to facilitate its removal after use. It is equipped with two adhesives (3) to ensure its maintenance on the patient's torso as well as a tracheal barrette (4). By pre-forming the tracheal barrette (4), a space is created between the tracheostoma and the device, thus preventing the protection from sticking to the tracheotomy and disturbing the patient's breathing.

COMPOSITION - The TRACHEOCLEAN® disposable tracheal protection is made of a non-woven material composed of 70% viscose and 30% polyester.

PLACEMENT BY A THIRD PARTY (NURSE, CARETAKER)

STEP 1: PLACEMENT OF THE DISPOSABLE TRACHEAL PROTECTION IN FRONT OF THE PATIENT'S TRACHEOSTOMA

Carefully wash your hands beforehand, then take the TRACHEOCLEAN® disposable tracheal protection by its two cords and place it in front of the patient's tracheostoma. The tracheal barrette and the two adhesives must be facing the patient.

STEP 2: ATTACHING THE DISPOSABLE TRACHEAL PROTECTION

Attach the TRACHEOCLEAN® disposable tracheal protection behind the patient's neck using the two cords and make a double knot.

STEP 3: PRE-FORMING OF THE TRACHEAL BARRETTE

Pre-form the tracheal barrette so that it creates a space between the TRACHEOCLEAN® disposable tracheal protection and the tracheal orifice.

STEP 4: PREPARATION OF THE ADHESIVES

In order to ensure that the TRACHEOCLEAN® disposable tracheal protection is maintained in place, particularly during the night, peel off the protection of the two adhesives.

STEP 5: PLACEMENT OF THE ADHESIVES

Stick the two adhesives on the torso by pressing gently

PLACEMENT BY THE PATIENT

REMOVAL OF THE TRACHEOCLEAN® DISPOSABLE TRACHEAL PROTECTION BY THE PATIENT OR BY A THIRD PARTY

If there is mucus on the TRACHEOCLEAN® disposable tracheal protection, remove the protection with the help of the pre-cut. To do this, pull on either side of the cord where the pre-cut has been made.

STOCKAGE

Store the TRACHEOCLEAN® in a dark, clean and dry place.

ELIMINATION

The TRACHEOCLEAN® disposable tracheal protection is a single-use device and must be disposed of either with treatment waste (hospital setting) or in the household waste bin (at the patient's home).

TRACHEOCLEAN® DISPOSABLE TRACHEAL PROTECTION				
Name	Ref.	LPPR code	Reimbursement by national health insurance system	Name to be entered
TRACHEOCLEAN (x100)	TRCLN-100	6184886	No limit	Tracheal protection for laryngectomy patients with stand-up collar (minimum order: 100 units)
TRACHEOCLEAN (x50)	TRCLN-050	6184886		

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